


## UPDATE BASED ON SOFTWARE VERSION 3.1.7

### 6.2.3 Graph


If you select the "Graph info" option, a graphical representation of the last 30 seconds of the mask pressure and the patient flow appears on the display in real time.



The current measured value is displayed on the left-hand side of each of the graphs.

The values of the therapy status area are displayed at the bottom (see Page 28).

In the Bilevel ST or APCV mode, the inspiration is indicated with S or T. (S) stands for spontaneous inspiration (T) stands for timed mode.



If any AHI events detected by the device, the corresponding event is indicated with A or H. (A) stands for apnea (H) stands for hypopnea.

## 6.3 Menu

### 6.3.1 Quick View

BACK	
PARAMETER	VALUE
DAILY THERAPY TIME	0 h/day
THERAPY HOURS	0 h
OPERATING TIME	17 min
AHI	0.0
SOFTWARE	3.1.7
FIRMWARE	3.0.15

MENU/QUICK VIEW  PSV / BILEVEL-ST  12:58

You will find a list of the therapy time used per day, the total usage of the device, the total operating time, the AHI of the current day and the current software/firmware version under Quick view.



#### **Warning**

*The AHI values shown in the Quick View is for reference only, cannot be used for diagnosis purpose.*

## 7.2 Handling errors

### 7.2.1 Error messages (display)

Error message	Possible cause	Measure
"The set temperature of the humidifier cannot be reached."	If the water filled in the water chamber has a significantly lower temperature than the room air, the heating time may last longer.	Fill in the water with room temperature. If the message persists, contact your service partner.
"Humidifier temperature is too high. Humidifier is deactivated."	<ul style="list-style-type: none"><li>- The heating mattress is defective.</li><li>- There is no water in the water chamber.</li></ul>	Unplug the device for 5 minutes. Restart the humidifier and check the water level in advance. If the message appears again, contact your service partner.
"The start temperature of the humidifier does not change."	<ul style="list-style-type: none"><li>- The heating mattress is defective.</li><li>- A non-compatible humidifier is connected.</li></ul>	Check if a compatible humidifier is connected. Contact your service partner if necessary.